



[Home](#) > [Stem Cell Lines](#) > [Order Stem Cell Lines](#) > **WA22**

WA22

Human embryonic stem (ES) cell line derived in MEF conditioned medium on Matrigel

Alias: N/A Cell Type: Human ES Culture Platform / Protocol: Feeder Independent - TeSR1 Medium Disease Model: N/A Genetic Modification Keyword: N/A NIH Registry Approved: Yes Karyotype: 46,XX Blood Type: O+ Publication: Publication not available Research Use Only: Yes	QTY 1-2 3-4 5+ PRICE \$1,250.00/val \$1,125.00/val \$1,062.50/val Quantity 1 ADD TO CART	<input type="checkbox"/> Check if you received this cell line previously 302x2
---	--	---

Current Lot Information

Lot Number	Lot Description	Passage Number	Banked By	Product Information and Testing
WR0056		12	WiCell	PDF

Historical Lot Information

Lot Number	Lot Description	Banked By	Product Information and Testing
WB0053		WiCell	PDF

EXPAND YOUR LAB WITH OURS

ORDER Stem Cell Lines Cytoprotective Services Stem Cell Services	Stem Cell Lines Order Stem Cell Lines About the WCell Stem Cell Bank Research and Clinical Applications Newly Available Adult Cell Lines Volume and Marketing Discounts Clinical Grade (cGMP) Cell Processes	Cytoprotective Services Send Samples for Testing Request Pricing Information Research and Clinical Applications Specific Handling (SKY) CDM/Manufacturing Fluorescence in Situ Hybridization (FISH) Short Tandem Repeat (STR) CDM+ Stem Monitoring Stem Microarray bioPDI	Stem Cell Services Banking Stem Cells Quality Control Testing Request Additional Information WCell Stem Cell Bank Technical Support Training and Outreach Webinars New Website Features	Support Order, Report and FAQs WCell Stem Cell Bank Technical Support Training and Outreach Webinars New Website Features	How WCell Plays Pivotal Contact Us Leadership Community Subscribe Privacy and Terms Newsroom
		Cytoprotective Reports and Privacy Policy Personalized Publications			

[Home](#) | [Newsroom](#) | [Contact Us](#) | [WABC](#) | [10W-Madison](#) | [Sitemap](#) | [Email a Friend](#) | [Print This Page](#)

WiCell Research Institute, Inc. is a nonprofit organization offering research and clinical grade pluripotent stem cell lines, cytogenetic testing, quality control testing and cell banking services to researchers worldwide.
Let WiCell be Your Lab Partner! ©2012 WiCell. All rights reserved.



Histocompatibility/Molecular Diagnostics Laboratory


Product Information and Testing - Amended

General Cell Line Testing Performed by WiCell

The following tests were performed on the cell line. The tests do not apply to any particular lot.

Test Description	Test Provider	Test Method
Differentiation Potential by Teratoma	WiCell	SOP-CH-213 SOP-CH-214
HLA	UW Molecular Diagnostics Laboratory	PowerPlex 1.2 System by Promega
ABO	American Red Cross	For ABO: Olsson ML, Chester MA. A rapid and simple ABO genotype screening method using a novel B/O2 versus A/O2 discriminating nucleotide substitution at the ABO locus. Vox Sang 1995; 69(3):242-7. For RHD: Singleton BK, Green CA, Avent ND, Martin PG, Smart E, Daka A, Narter-Olaga EG, Hawthorne LM, Daniels G. The presence of an RHD pseudogene containing a 37 base pair duplication and a nonsense mutation in Africans with the Rh D-negative blood group phenotype. Blood 2000; 95(1): 12-8.
Growth Curve (Doubling Time)	WiCell	Varies by culture platform
Flow Cytometry for ESC Marker Expression	UW Flow Cytometry Laboratory	SOP-CH-101 SOP-CH-102 SOP-CH-103 SOP-CH-105
Array Comparative Genomic Hybridization (aCGH)	WiCell	SOP-CH-308 SOP-CH-309 SOP-CH-310
Comprehensive Human Virus Panel	Charles River	ID910

Amendment(s)	Reason for Amendment	Date
CoA updated to include copyright information.		See Signature
CoA updated for format changes, including adding fields of thaw recommendation, vial label, protocol, and banked by, and removal of footnotes. General Cell Line Testing CoA added to lot CoA.		24-JUN-2013
Original CoA		18-MAR-2011

Date of Lot Release	Quality Assurance Approval
18-March-2011	<div style="text-align: right;">1/3/2014</div> <div style="text-align: center;">  AMC Quality Assurance </div>

Short Tandem Repeat Analysis*

Sample Report: 8155-STR

UW HLA#: 63885

Sample Date: 10/01/10
Received Date: 10/01/10

Requestor: WiCell Research Institute

Test Date: 10/05/10

File Name: 101005SLE

Report Date: 10/06/10

Sample Name: (label on tube) 8155-STR

Description: WiCell Research Institute provided genomic DNA
248.5 µg/mL; 260/280 = 1.93

Locus	Repeat #	STR Genotype
D16S539	5, 8-15	11,14
D7S820	6-14	10,11
D13S317	7-15	12,12
D5S818	7-15	13,13
CSF1PO	6-15	11,12
TPOX	6-13	8,9
Amelogenin	NA	X,X
TH01	5-11	6,6
vWA	11, 13,21	17,19

Comments: Based on the 8155-STR DNA data and received on 10/01/10 from WiCell Research Institute, this sample (U/LH/AF-63885) exactly matches the STR profile of the human stem cell line W422 comprising 13 allelic polymorphisms across the 8 STR loci analyzed. No STR polymorphisms other than those corresponding to the human W422 stem cell line were detected and the concentration of DNA required to achieve an acceptable STR genotype (signal/noise) was equivalent to that required for the standard procedure (~1 ng/amplification reaction) from human genomic DNA. These results suggest that the 8155-STR DNA sample submitted corresponds to the W422 stem cell line and it was not contaminated with any other human stem cells or significant non-human animal contamination. Sensitivity limits for detection of STR polymorphisms unique to either this or other human stem cell lines is ~5%.

HLA/Molecular Diagnostics Laboratory

HLA/Molecular Diagnostics Laboratory

Report Number
852056.A01
Page 1 of 1December 1, 2010
P.O. # [REDACTED]
AMENDED REPORT
Original Issue Date:
11-27-10

Amendment Summary

WuXi Research Institute

STERILITY TEST REPORT

Sample Information:

HES Cells

- 1: WA15-07-07-WB0062 #1661
- 2: WA22-WB0046 #1491
- 3: WA13-C-WB0054 #7289
- 4: WA22-WB0053 #3855
- 5: IPS(WR50)-3-WB0057 #3060
- 6: WA23-WB0067 #4096
- 7: WA15-07-03-WB0063 #8295

Date Received:
Date in Test:
Date Completed:November 09, 2010
November 11, 2010
November 25, 2010

Test Information:

Test Codes: 30744, 30744A
Immersion, USP / 21 CFR 610.12
Procedure #: BS210WCR.201

TEST PARAMETERS	PRODUCT	
Approximate Volume Tested	0.5 mL	0.5 mL
Number Tested	14	14
Type of Media	SCD	FTM
Media Volume	400 mL	400 mL
Incubation Period	14 Days	14 Days
Incubation Temperature	20 °C to 26 °C	30 °C to 35 °C
RESULTS	14 NEGATIVE	14 NEGATIVE

AD1 - Dated 12-01-10: Corrected sample information for sample #1.

QA Reviewer
Date

12-01-10

Technical Reviewer
Date

12-01-10

Testing conducted in accordance with current Good Manufacturing Practices.

APPENDIX

Document ID #: DCF9002E
Title: QUALITY ASSURANCE REPORT - GMP
Effective Date: 03/24/10
Edition #: 03

QUALITY ASSURANCE REPORT - GMP

TEST PERFORMED	PROCEDURAL REFERENCE	TEST PERFORMED	PROCEDURAL REFERENCE
<input checked="" type="checkbox"/> M-250	SOP's 3008, 3011, 3013	<input type="checkbox"/> M-700	SOP's 3008, 3009, 3010
<input type="checkbox"/> M-300	SOP's 3008, 3014	<input type="checkbox"/> M-800	SOP's 3008, 3011, 3016
<input type="checkbox"/> M-350	SOP's 3008, 3014, 3015		

Bionique Sample ID #(s) 62621

This testing procedure was performed in compliance with the FDA's Current Good Manufacturing Practice (cGMP) standards (to the extent that the regulations pertain to the procedures performed) as specified in the Code of Federal Regulations, Title 21 Parts 210 and 211 [21 CFR 210 & 211]. All related records derived from the test procedures have been reviewed by the Quality Assurance Department. The individual's signature below verifies that the methods and procedures referenced above have been followed and that the Final Report accurately reflects the raw data generated during the course of the procedures. All records, including raw data and final reports are archived on site for a minimum of seven years.

The specified test's procedures determine the intervals at which samples are inspected. The medium used for testing must pass quality control mycoplasma growth promotion testing and sterility testing. Traceability of all of the components used is assured and supporting documentation can be supplied upon request.

Quality Assurance Review Date: 10/27/10

Reviewed By: [REDACTED] A Associate [REDACTED]

NOTE:

1. Prior to receipt at Bionique® Testing Laboratories, Inc., the stability of the test article is the responsibility of the company submitting the sample. Bionique Testing Laboratories Inc. will assume responsibility for sample stability following receipt and prior to being placed on test.
2. This test is for the detection of microbiological growth and does not require statistical validation.

Page 1 of 2

APPENDIX

BIONIQUE® TESTING LABORATORIES, INC.

Document ID #: DCF9002E
Title: QUALITY ASSURANCE REPORT - GMP
Effective Date: 03/24/10
Edition #: 03

REFERENCES

Regulatory:

1. Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Code of Federal Regulations [CFR], Title 21 CFR Part 210, Current Good Manufacturing Practice in Manufacturing, Processing, Packing, or Holding of Drugs; General. FDA. Office of the Federal Register, National Archives and Records Department.
2. Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Code of Federal Regulations [CFR], Title 21 CFR Part 211, Current Good Manufacturing Practice for Finished Pharmaceuticals. FDA. Office of the Federal Register, National Archives and Records Department.
3. Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Points to Consider in the Characterization of Cell Lines Used to Produce Biologicals, Director, Center for Biologics Evaluation and Research, May, 1993. Docket No. 84N-0154.
4. Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Code of Federal Regulations [CFR], Title 21 CFR Part 610.30, General Biological Products Standards; Subpart D, Test for Mycoplasma. FDA. Office of the Federal Register, National Archives and Records Department.

General:

1. Barile MF, Kern J. Isolation of Mycoplasma arginini from commercial bovine sera and its implication in contaminated cell cultures. Proceedings of the Society for Experimental Biology and Medicine, Volume 138, Number 2, November 1971.
2. Chen, T.R. In situ detection of mycoplasma contamination in cell cultures by fluorescent Hoechst 33258 stain. Experimental Cell Research, 104: 255-262, 1977.
3. Carolyn K. Lincoln and Daniel J. Lundin. Mycoplasma Detection and Control. U. S. Fed. for Culture Collections Newsletter, Vol. 20, Number 4, 1990.
4. Fetal Bovine Serum; Proposed Guideline. National Committee For Clinical Laboratory Standards (NCCLS), Vol. 10, Number 6, 1990. (NCCLS publication M25-P).
5. McGarrity GJ, Sarama J, Vanaman V. Cell Culture Techniques. ASM News, Vol. 51, No. 4, 1985.
6. Tully JG, Razin S. Methods in Mycoplasmaology, Volumes I and II. Academic Press, N.Y., 1983.
7. Barile MF, Razin S, Tully JG, Whitcomb RF. The Mycoplasmas, Volumes 1-4. Academic Press, N.Y., 1979.
8. <http://www.bionique.com/> - Safe Cells Insights

APPENDIX IV

Page 1 of 2

Document#: DCF3013D
Edition#: 10
Effective Date: 07/15/2003
Title: M-250 FINAL REPORT SHEET

M-250 FINAL REPORT

Direct Spindown Culture
Procedure 3008, 3011, 3013TO: WuCell QA
WuCell Research Institute

BTL SAMPLE ID#: 62621 P.O.#: [REDACTED] DATE REC'D: 09/28/2010

TEST/CONTROL ARTICLE:

WA22-WB0053 #8155

LOT#: NA

DIRECT CULTURE SET-UP (DAY 0) DATE: 09/29/2010
INDICATOR CELL LINE (VERO) SEE DIA FLOUROCHROME MICRO SHEET

			DATE
THIOGLYCOLLATE BROTH	DAY 7	+	10/06/2010
	DAY 28	+	10/27/2010
BROTH-FORTIFIED COMMERCIAL			
0.5 mL SAMPLE	DAY 7	+	10/06/2010
6.0 mL BROTH	DAY 28	+	10/27/2010
BROTH-MODIFIED HAYFLICK			
0.5 mL SAMPLE	DAY 7	+	10/06/2010
6.0 mL BROTH	DAY 28	+	10/27/2010
BROTH-HEART INFUSION			
0.5 mL SAMPLE	DAY 7	+	10/06/2010
6.0 mL BROTH	DAY 28	+	10/27/2010

(See Reverse)

Document#: DCF3013D
Edition#: 10
Effective Date: 07/15/2003
Title: M-250 FINAL REPORT SHEET

SAMPLE ID#:	62621	AEROBIC	MICROAEROBIC	DATE
AGAR PLATES-FORTIFIED COMMERCIAL	DAY 7 DAY 14 DAY 21	+ + +	+ + +	10/06/2010 10/13/2010 10/20/2010
AGAR PLATES-MODIFIED HAYFLICK	DAY 7 DAY 14 DAY 21	+ + +	+ + +	10/06/2010 10/13/2010 10/20/2010
AGAR PLATES-HEART INFUSION	DAY 7 DAY 14 DAY 21	+ + +	+ + +	10/06/2010 10/13/2010 10/20/2010
BROTH SUBCULTURES (DAY 7)				
	DATE: 10/06/2010			
AGAR PLATES-FORTIFIED COMMERCIAL	DAY 7 DAY 14 DAY 21	+ + +	+ + +	10/13/2010 10/20/2010 10/27/2010
AGAR PLATES-MODIFIED HAYFLICK	DAY 7 DAY 14 DAY 21	+ + +	+ + +	10/13/2010 10/20/2010 10/27/2010
AGAR PLATES-HEART INFUSION	DAY 7 DAY 14 DAY 21	+ + +	+ + +	10/13/2010 10/20/2010 10/27/2010

RESULTS: No detectable mycoplasma contamination

Date: 10/27/10
Laboratory Director: [Redacted] Ph.D.

W-250 Procedure Summary: The objective of this test is to determine whether or not detectable mycoplasmas are present in an in vitro cell culture system. In a primary culture, mycoplasma would show as cell lines. This procedure involves an indirect, non-invasive approach to detect non-culturable mycoplasmas with a direct culture methodology utilizing three different mycoplasma media formulations. The indirect approach involves the inoculation of the sample into a mycoplasma-free VERO (Vero) indicator cell line and performing a DNA fluorescence assay after 10-120 hours of incubation. The direct culture aspect of the test utilizes three different mycoplasma media including both direct and agar formulations. The ability to transcribe DNA from the 3' end of the 16S rRNA gene and the ability to amplify the 16S rRNA gene of the 3' end of the 16S rRNA gene. Samples from both to direct and direct to agar media are used. After incubation, agar media are observed visually and microscopically in order to detect any visible growth and morphologically indicative of mycoplasma contamination. Absence of the final result with respect to the laboratory director's report that the sample material was confirmed consistently with the test procedure as detailed in the referenced SOP and that all test conditions have been found to meet the required acceptance criteria for a valid test, including the appropriate results for the positive and negative controls.

Document ID #: DCF3008A
Title: DNA FLUOROCHROME ASSAY RESULTS
Effective Date: 3/24/10
Edition #: 07

DNA-FLUOROCHROME ASSAY RESULTS
Procedures 3008, 3009, 3011

Sample ID #: 62621 M-250 Date Rec'd: 09/28/2010 P.O. #: [Redacted]
Indicator Cells Inoculated: Date/Initials: 9/30/10 / HB
Fixation: Date/Initials: 10/4/10 / HB
Staining: Date/Initials: 10/4/10 / HB

TEST/CONTROL ARTICLE:

WA22-WB0053 #8155

LOT# NA

WiCell QA
WiCell Research Institute

Phone: [Redacted]

Fax #: [Redacted]

DNA FLUOROCHROME ASSAY RESULTS:

☒ **NEGATIVE:** A reaction with staining limited to the nuclear region, which indicates no mycoplasma contamination.

☐ **POSITIVE:** A significant amount of extranuclear staining which strongly suggests mycoplasma contamination.

☐ **INCONCLUSIVE:**
A significant amount of extranuclear staining consistent with low - level mycoplasma contamination or nuclear degeneration.
A significant amount of extranuclear staining consistent with bacterial, fungal or other microbial contaminant or viral CPE. Morphology not consistent for mycoplasma contamination.

COMMENTS:

Date: 10/4/10 Results Read by: HB Date of Review: 10/4/10 Reviewed by: Self



WiCell Cytoogenetics Report: 003688
WISC 8155

Report Date: September 26, 2010

Case Details:

Cell Line: WA22-WB0053 (8155)

Passage #: 11

Date Completed: 9/26/2010

Cell Line Gender: Female

Investigator: Wisconsin International Stem Cell Bank

Specimen: hESC on Matrigel

Date of Sample: 9/20/2010


Tests/Reason for: lot release testing

Results: 46,XX

Completed by: [Redacted] CG(ASCP), on 9/23/2010

Reviewed and interpreted by: [Redacted] PhD, FACMG, on 9/26/2010

Interpretation: No abnormalities were detected at the stated band level of resolution.



Cell: S01-03
Slide: 2-13
Slide Type: Karyotyping

of Cells Counted: 20
of Cells Karyotyped: 4
of Cells Analyzed: 8
Band Level: 450-500

Results Transmitted by Fax / Email / Post
Sent By: [Redacted]
QC Review By: [Redacted]

Date: [Redacted]
Sent To: [Redacted]
Results Recorded: [Redacted]

	Teratoma Histology Report	FORM SOP-CH-014.01 Version B Edition 01
--	---------------------------	---

Cell Line: WA22 Cell Lot Number: NA Sample Number: 9971

ECTODERM	
Structure Name: Brain Magnification: 200X Slide ID: A	Structure Name: Neuroendocrine Magnification: 200X Slide ID: A
ENDODERM	
Structure Name: Hepatoid Magnification: 200X Slide ID: A	Structure Name: Bronchial mucosa Magnification: 200X Slide ID: A
MESODERM	
Structure Name: Cartilage Magnification: 200X Slide ID: A	Structure Name: Nephroid Magnification: 200X Slide ID: B

Comments: Structures identified include Ectoderm (2), Mesoderm (2) and Endoderm (2)

Sample(s) were assessed for the presence of differentiation into cell types characteristic of the three embryonic germ layers, which, if present in the sample(s) examined, are represented in the photographs above. The individual's signature below verifies that this report accurately reflects the pathology observed.

Pathologist (By/Date): [Redacted]

QA Review (By/Date): [Redacted]

Date: 09/02/2010 17:10:33

To: Wicell Research Institute

Re: High-resolution HLA results

Patient

Name		HLA DNA-based typing*								
HLA / MR# received	Dates	Method:	PCR-SSP	Direct Sequencing						PCR-SSP
		A*	B*	C*	DRB1*	DRB3*	DRB4*		DRB5*	DOB1*
WICELL 8432-HLA	DOB SNP	02.01	14.02	03.04	01.02					
63679	A,B,C SNP	09/02/2010	MR02	40.01	08.02	08.01				
09/02/2010	DRB Seq	09/02/2010								

December 9, 2010

Wicell Research Institute

SAMPLE: DNA WA22 8432 (MA#388-10)

Date Received: 11/17/10

Sample Date: 08/26/10

HISTORY: DNA from cell line.

TEST REQUESTED: Genotype for ABO and common RH

TESTING PERFORMED: ABO: Polymerase chain reaction-restriction fragment length polymorphism (PCR-RFLP) testing for nucleotide (nt) positions 261 (O⁺), 467 (A⁺), 703 (B), and 1096 (B and O⁺). RH: Multiplex PCR-RFLP for RHD and RHCE^{C/c}. HEA Beadchip for RHCE^{E/e}.

DNA RESULTS: PCR-RFLP indicated heterozygous for at 261G characteristic of O⁺ alleles.

Result	Test Method
ABO ⁺ O ⁺ O ⁺	PCR-RFLP
RHD positive for exons 4, 7 and no inactivating pseudogene	Multiplex PCR
RHCE ^{C/c}	Multiplex PCR
RHCE ^{E/e}	HEA 2.1 Assay

Predicted phenotype: Group O, RhD+C-E++

Manager, Molecular Analysis

Director, Immunohematology and Genomics

These *in vitro* diagnostic tests were developed and their performance characteristics established in the Molecular Analysis Laboratory. The tests have not been submitted to the Food and Drug Administration (FDA) for clearance or approval and, therefore, are not FDA-licensed tests. The Molecular Analysis Laboratory is certified under the Clinical Laboratory Improvement Amendments (CLIA) of 1988 as qualified to perform high complexity clinical testing. The New York Blood Center has been approved, by the New York State Department of Health to perform these tests under its current Clinical Laboratory Permit. These results are intended to predict a blood group antigen profile in a patient or donor, and are not intended for clinical diagnosis or as the sole means for patient management decisions. There are situations where testing DNA of a person may not reflect the red cell phenotype and not all performance characteristics have been determined. Nucleotide changes that inactivate gene expression or rare new variant alleles may not be identified in these assays.



HLA/Molecular Diagnostics Laboratory

Date:



HLA/Molecular Diagnostics Laboratory

Date:


This test was developed and its performance characteristics determined by the UWMC Clinical Laboratory. It has not been cleared or approved by the U.S. Food and Drug Administration. However, the FDA does not regulate because of analytic specific reagents since the laboratory is approved, under CLIA, for high complexity testing.

©2009 Wisconsin Hospital and Clinics



Characterization Report- Growth Characteristics

Cell Line Information	NSCB QA Use
Sample ID: 4122	Cell lot #: New Derivation
Cell Line: WA22-A in mTeSR1	Report reviewed by: JKT
Passage: p12	Report reviewed on: 13Oct10
	Date cells received: 17Aug10



Doubling time and confidence Interval data:

Slope ± 95% C.I. 0.03596 ± 0.002378

Doubling Time ± 95% C.I.


20.55 hours ± 1.4 hours=

19.15 hours ~ 21.95 hours



Characterization Report- Growth Characteristics



Cell Line Information	NSCB QA Use
Sample ID: 4122	Cell lot #: New Derivation
Cell Line: WA22-A in mTeSR1	Report reviewed by: JKT
Passage: p12	Report reviewed on: 13Oct10
	Date cells received: 17Aug10





Procedures performed: SOP-CH-101 SOP-CH-102 SOP-CH-103 SOP-CH-105
Cell Line: WA22 TeSR/MG Passage 13 Sample ID: 4122-FAC
Date of: (mm/dd/yy) acquisition: 09/17/10 file creation: 09/17/10 file submission: 09/20/10

antigen2:	SSEA4 - antigen2 +	SSEA4 + antigen2 +	SSEA4 + antigen2 -	SSEA4 - antigen2 -	ALL SSEA4 +	ALL antigen2 +
SSEA3	0.33	99.10	0.38	0.22	99.48	99.43
TRA1-60	0.69	98.20	1.11	0.02	99.31	98.89
TRA1-81	0.20	99.00	0.81	0.00	99.81	99.20
Oct-4	0.52	96.30	2.94	0.21	99.24	96.82
SSEA1	0.00	0.79	99.10	0.06	99.89	0.79



WiCell CytoGenetics Report: 003689 WISC2100

Report Date: 7/1/2011
Date of Sample: 9/24/2010
Investigator: [REDACTED]
Reason for Testing: lot release testing
Specimen: hESC on Matrigel, TeSR
Karyotype Results: n/a

Test: WA22-WB0046p10 (Female)
Reference: WA01-MCB-03-S-5p26(3) (Male)
Project: 221
Funding: 000
CGH Accession #: 000398
GEO Accession #: [REDACTED]

Microarray Results:
☒ arr(1-22,X)2 - Female
☐ arr(1-22)x2,(XY)x1 - Male
☐ Consistent with the Karyotype Results
☐ Inconsistent with the Karyotype Results
☐ Consistent with a Balanced Karyotype (Karyotype Unavailable)
☐ Additional Findings

Interpretation:
CNV gains/losses
There were 34 copy number gains and losses identified, including 2 pseudoautosomal regions and 8 copy number changes due to the reference DNA
Select CNVs are detailed in the table below

Chr	Band (Genomic Position)	Width	Aberration Type	Classification	Genes
1	arr 1q42.3(232,994,064-233,017,150)x1	23,086	Loss	Uncertain Significance - Likely Benign	
1	arr 1q43(241,139,770-241,196,609)x1	56,838	Loss	Uncertain Significance - Likely Benign	
2	arr 2q37.3(242,535,552-242,648,925)x1	113,372	Loss	Uncertain Significance - Likely Benign	
7	arr 7p13(43,966,453-44,047,927)x1	81,474	Loss	Uncertain Significance - Likely Benign	DBNL, UBE2D, WSCR19
7	arr 7q22.1(101,904,922-102,096,488)x1	191,566	Loss	Uncertain Significance - Likely Benign	LRWD1, MGC119295, POLR21, POLR22, POLR2J3, BASA4
7	arr 7q35(143,306,579-143,705,123)x3	398,544	Gain	Uncertain Significance - Likely Benign	ABHGE5, FLH3692, OR2A1, OR2A12, OR2A14, OR2A2, OR2A25, OR2A42, OR2A5, OR2A1, OR6B1
9	arr 9p23(12,111,305-12,361,968)x1	250,662	Loss	Uncertain Significance - Likely Benign	
10	arr 10q26.3(135,102,844-135,187,332)x3	84,487	Gain	Uncertain Significance - Likely Benign	CYP2E1
12	arr 12q24.1(113,781,059-113,814,033)x1	32,974	Loss	Uncertain Significance - Likely Benign	
17	arr 17p11.2(18,303,144-18,349,050)x1	45,906	Loss	Uncertain Significance - Likely Benign	LOC654346
17	arr 17q21.31q21.32(41,709,705-42,238,590)x1	528,884	Loss	Uncertain Significance - Likely Benign	ARL17, ARL17P1, LRR37A, LRR37A2, NSF, WNT3
19	arr 19q13.33(56,832,782-56,853,080)x1	20,298	Loss	Uncertain Significance - Likely Benign	SIGLEC14, SIGLEC5
19	arr 19q13.42(59,231,079-59,251,060)x1	19,981	Loss	Uncertain Significance - Likely Benign	VSTM1

Notes:
Karyotype Information - n/a
Published CNVs (4) - Narva et al: arr 15q11.2(18,469,957-20,226,623)x3

References: Werbowitski-Ogilvie, T., Bosse, M., Stewart, M., et al. (2008). Characterization of human embryonic stem cells with features of neoplastic progression. *Nature Biotechnology* 27, 91-97.
Wu, H., Kim, K., Mehra, K., et al. (2008). Copy number variant analysis of human embryonic stem cells. *Stem Cells* 26, 1484-1489.
Chin, M.H., Mason, M., Xie, W., et al. (2009). Induced pluripotent stem cells and embryonic stem cells are distinguished by gene expression signatures. *Cell Stem Cell* 5, 111-123.
Narva, E., Auto, R., Rahkonen, N., et al. (2010). High-resolution DNA analysis of human embryonic stem cell lines reveals culture-induced copy number changes and loss of heterozygosity. *Nature Biotechnology* 28, 371-377

Page | 1

Recommendations: For relevant findings, confirmation and localization is recommended. Contact cytogenetics@wicell.org to request further testing.

Results Completed By: [REDACTED]
Reviewed and Interpreted By: [REDACTED]

aCGH Specifications:

- Platform: NimbleGen 12x135K array (HG18 WG CGH v3.1 HX12)
- Relative copy number is determined by competitive differential hybridization of labeled genomic DNA to the 135,000 oligonucleotide whole genome tiling array
- Probe length = 60mer, spanning non-repetitive regions of the human genome
- Median probe spacing = 21,500
- Analysis software: NimbleScan™, CGH Fusion (RBS v1.0)™
- Array design, genomic position, genes and chromosome banding are based on HG18
- Analysis is based on examination of unaveraged and/or 130Kbp (10K) averaged data tracks as noted. Settings for data analysis in Infoquant include an average log-ratio threshold of 0.2, a minimum aberration length of 5 probes, p-value of 0.001. Additional analysis of this data may be performed using different ratio settings and different window averaging to enhance resolution.
- Raw data has not yet been deposited in GEO.
- Reported gains and losses are based on test to reference ratios within CGHfusion™ and the size of aberration.
- Quality assurance monitors: 1) opposite gender reference DNA ratio change in X and Y chromosomes; 2) presence of Xpter and Xq21.3 "pseudoautosomal" (PAR) imbalance; 3) presence of known reference DNA copy number changes. QA measures—PAR (2/2). Reference DNA copy number changes (B); test sample gain or loss of X and Y chromosomes consistent with the opposite gender reference sample.

Limitations: This assay will detect aneuploidy, deletions, duplications of represented loci, but will not detect balanced alterations (reciprocal translocations, Robertsonian translocations, inversions, and insertions), point mutations, loss of heterozygosity (LOH), uniparental disomy or imbalances less than 30kb in size. Copy number variants can be attributable to the test or reference samples used. Exact limits of detectable mosaicism have not been determined, but >20% mosaicism is reported to be visualized by aCGH. Actual chromosomal localization of copy number change is not determined by this assay. Other mapping procedures are required for determining chromosomal localization.

Results Transmitted by ☐ Fax / ☐ Email / ☐ Post
Sent By: [REDACTED]
Date: [REDACTED]
Sent To: [REDACTED]

Printed: Thursday, November 18, 2010 at 9:3

Charles River Research Animal Diagnostic Services

Sponsor: WiCell Research Institute

Accession #: 2010-048114

Diagnostic Summary Report

Received: 16 Nov 2010
Approved: 18 Nov 2010, 09:30
Bill Method: PO# [REDACTED]
Test Specimen: Human

Sample Set	Service (# Tested)	Profile	Assay	Tested	+	+/-	?
#1	Infectious Disease PCR (3)	All Results Negative					

+ = Positive, +/- = Equivocal, ? = Indeterminate

Service Approvals

Service	Approved By*	Date
Infectious Disease PCR	[REDACTED]	18 Nov 2010, 09:27

To assure the SPF status of your research animal colonies, it is essential that you understand the sources, pathobiology, diagnosis and control of pathogens and other adventitious infectious agents that may cause research interference. We have summarized this important information in infectious agent **Technical Sheets**, which you can view by visiting http://www.criver.com/info/disease_sheets

*This report has been electronically signed by laboratory personnel. The name of the individual who approved these results appears in the header of this service report. All services are performed in accordance with and subject to General Terms and Conditions of Sale found in the Charles River Laboratories-Research Models and Services catalogue and on the back of invoices.

CR-RADS-ILIMS Form: FM-1741 Rev: 3

Page 1 of 2

Sponsor: WiCell Research Institute

Accession #: 2010-048114

Product: Not Indicated

Test Specimen: Human

Received: 16 Nov 2010

Molecular Diagnostics Infectious Disease PCR Results Report

Department Review: Approved by [Signature] 18 Nov 2010, 09:27*

Human Comprehensive Virus Panel

Sample #: Code :	1 WA22-WB0046 49128	2 WA22-WB0047 49010	3 WA22-WB0046 49032
John Cunningham virus	-	-	-
BK virus	-	-	-
Herpesvirus type 6	-	-	-
Herpesvirus type 7	-	-	-
Herpesvirus type 8	-	-	-
Parvovirus B19	-	-	-
Epstein-Barr Virus	-	-	-
Hepatitis A virus	-	-	-
Hepatitis B virus	-	-	-
Hepatitis C virus	-	-	-
HPV-16	-	-	-
HPV-18	-	-	-
Human T-lymphotropic virus	-	-	-
Human cytomegalovirus	-	-	-
HIV-1	-	-	-
HIV-2	-	-	-
Adeno-associated virus	-	-	-
Human Foamy Virus	-	-	-
LCMV PCR	-	-	-
Hantavirus Hantaan PCR	-	-	-
Hantavirus Seoul PCR	-	-	-
Mycoplasma Genus PCR	-	-	-
DNA Spike	PASS	PASS	PASS
RNA Spike	PASS	PASS	PASS
NRC	PASS	PASS	PASS

Remarks: - = Negative; I = Inhibition; +/- = Equivocal; + = Positive.

Sample Suitability/Detection of PCR Inhibition:

Sample DNA or RNA is spiked with a low-copy number of an exogenous DNA or RNA template respectively. A spike template-specific PCR assay is used to test for the spike template for the purpose of determining the presence of PCR inhibitors. The RNA spike control is also used to evaluate the reverse-transcription of RNA. Amplification of spike template indicates that there is no detectable inhibition and the assay is valid.

NRC:

The nucleic acid recovery control (NRC) is used to evaluate the recovery of DNA/RNA from the nucleic acid isolation process. The test article is spiked with a low-copy number of DNA/RNA template prior to nucleic acid isolation. A template-specific PCR assay is used to detect the DNA/RNA spike.

*This report has been electronically signed by laboratory personnel. The name of the individual who approved these results appears in the header of this service report.

CR RADSLIMS Form: FM-1741 Rev. 3

Page 2 of 2



Product Information and Testing Amended Product Information

Product Name	WA22
Lot Number	WB0056
Parent Material	WA22-WB0046
Depositor	WiCell
Banked by	WiCell
Thaw Recommendation	Thaw 1 vial into 4 wells of a 6 well plate.
Culture Platform	Feeder Independent Medium: mTeSR1 Matrix: Matrigel
Protocol	WiCell Feeder Independent Protocol
Passage Number	p12 These cells were cultured for 11 passages prior to freeze. Cells were derived in Conditioned Medium on Matrigel. They were transitioned to mTeSR1 at passage 6 and cultured 5 additional passages prior to freeze. WiCell adds +1 to the passage number at freeze so that the number on the vial best represents the overall passage number of the cells at thaw.
Date Vial	29 September 2010
Vial Label	WB0056 WA22 p12 MW 29SEPT10
Biosafety and Use Information	Appropriate biosafety precautions should be followed when working with these cells. The end user is responsible for ensuring that the cells are handled and stored in an appropriate manner. WiCell is not responsible for damages or injuries that may result from the use of these cells. Cells distributed by WiCell are intended for research purposes only and are not intended for use in humans.

Testing Performed by WiCell

Test Description	Test Provider	Test Method	Test Specification	Result
Post-Thaw Viable Cell Recovery	WiCell	SOP-CH-305	≥ 15 Undifferentiated Colonies, ≤ 30% Differentiation	Pass
Identity by STR	UW Molecular Diagnostics Laboratory	PowerPlex 1.2 System by Promega	Consistent with known profile	Pass
Sterility - Direct transfer method	Appligene	30744	Negative	Pass
Mycoplasma	Bionique	M250	No contamination detected	Pass
Karyotype by G-banding	WiCell	SOP-CH-003	Normal karyotype	Pass

Amendment(s):

Reason for Amendment	Date
CoA updated for format changes, inclusion of banked by, protocol, and vial label.	See signature
Original CoA	05-April-2011

Date of Lot Release	Quality Assurance Approval
05-April-2013	 X AMC Quality Assurance Date: [Signature]

The material provided under this certificate has been subjected to the tests specified and the results and data described herein are accurate based on WiCell's reasonable knowledge and belief. Appropriate Biosafety Level practices and universal precautions should always be used with this material. For clarity, the foregoing is governed solely by WiCell's Terms and Conditions of Service, which can be found at <http://www.wicell.org/terms-conditions>.

© 2013 WiCell Research Institute



Histocompatibility/Molecular Diagnostics Laboratory
D4/231; (608) 263-8815
600 Highland Avenue
Madison, WI 53792-2472

Short Tandem Repeat Analysis®

Sample Report: 10027-STR

UW HLA#: 64698

Sample Date: 02/18/11

Received Date: 02/18/11

Requester: WiCell Research Institute

Test Date: 02/22/11

File Name: 110222 hb

Report Date: 02/23/11

Sample Name: (label on tube) 10027-STR

Description: WiCell Research Institute
provided genomic DNA
95.11 ug/mL; 280/280 = 1.93

Locus	Repeat #	STR Genotype
D10S1249	5, 6-13	11, 14
D7S820	6-14	10, 11
D10S1247	7-13	12, 12
D15S818	7-13	13, 13
CSF1PO	6-13	11, 12
TPX3X	6-13	8, 9
Adenogonin	NA	X, X
TH01	5-11	6, 6
+90A	11, 13, 21	17, 19

Comments: Based on the 10027-STR DNA dated and received on 02/18/11 from WiCell Research Institute, this sample (UW HLA# 64698) exactly matches the STR profile of the human stem cell line WA22 comprising 13 allelic polymorphisms across the 8 STR loci analyzed. No STR polymorphisms other than those corresponding to the human WA22 stem cell line were detected and the concentration of DNA required to achieve an acceptable STR genotype (signal/noise) was equivalent to that required for the standard procedure (1-1 ng/amplification reaction) from human genomic DNA. These results suggest that the 10027-STR DNA sample submitted corresponds to the WA22 stem cell line and it was not contaminated with any other human stem cells or a significant amount of mouse feeder layer cells. Sensitivity limits for detection of STR polymorphisms unique to either this or other human stem cell lines is ~5%.

Keith Chalotter, Manager
Molecular Diagnostics Laboratory

William M. Rebrauer, PhD, Director
Molecular Diagnostics Laboratory

Test Facility:
1265 Kennestone Circle
Marietta, GA 30066

This report is confidential. You may not be
allowed for reworking or other, unauthorized
without written permission. Results being sent
to the Laboratory Director.



WiCell Research Institute 505 S. Rosa Road Suite 120 Madison, WI 53716 Attn: Jessica Martin	Report Number 862282 Page 1 of 1 March 30, 2011 P.O. #: RP3934
---	--

STERILITY TEST REPORT

Sample Information:

HE5 Cells
1. WA22-WB0056 10059
2. WA21-WB0051 10060
3. WA24-WB0079 10061
4. WA07-WB0081 10062

Date Received:

March 10, 2011

Date in Test:

March 15, 2011

Date Completed:

March 29, 2011

Test Information:

Test Codes: 30744, 30744A
Immersion, USP 1 21 CFR 610.12
Procedure #: B5210WCR 201

TEST PARAMETERS	PRODUCT	
Approximate Volume Tested	0.5 mL	0.5 mL
Number Tested	8	8
Type of Media	SCD	FTM
Media Volume	400 mL	400 mL
Incubation Period	14 Days	14 Days
Incubation Temperature	20 °C to 25 °C	30 °C to 35 °C
RESULTS	8 NEGATIVE	8 NEGATIVE

GA Reviewer
Date

Date

Technical Reviewer
Date

Testing conducted in accordance with current Good Manufacturing Practices.



1265 Kennestone Circle • Marietta, GA 30066 • 888.847.9033 • 770.514.0262 • Fax 770.514.0294

* Testing to assess engraftment following bone marrow transplantation was accomplished by analysis of human genetic polymorphisms at STR loci. This methodology has not yet been approved by the FDA and is for investigational use only.

File: Final STR Report

APPENDIX

Document ID #: DCF9002F
Title: QUALITY ASSURANCE REPORT - GMP
Effective Date: 03/12/10
Edition #: 01

QUALITY ASSURANCE REPORT - GMP

TEST PERFORMED	PROCEDURAL REFERENCE	TEST PERFORMED	PROCEDURAL REFERENCE
<input checked="" type="checkbox"/> M-250	SOP's 3008, 3011, 3013	<input type="checkbox"/> M-700	SOP's 3008, 3009, 3010
<input checked="" type="checkbox"/> M-300	SOP's 3008, 3014	<input type="checkbox"/> M-800	SOP's 3008, 3011, 3016
<input type="checkbox"/> M-350	SOP's 3008, 3014, 3015		

Bionique Sample ID #(s) 64139

This testing procedure was performed in compliance with the FDA's Current Good Manufacturing Practice (cGMP) standards (to the extent that the regulations pertain to the procedures performed) as specified in the Code of Federal Regulations, Title 21 Parts 210 and 211 [21 CFR 210 & 211]. All related records derived from the test procedures have been reviewed by the Quality Assurance Department. The individual's signature below verifies that the methods and procedures referenced above have been followed and that the Final Report accurately reflects the raw data generated during the course of the procedures. All records, including raw data and final reports are archived on site for a minimum of seven years.

The specified test's procedures determine the intervals at which samples are inspected. The medium used for testing must pass quality control mycoplasma growth promotion testing and sterility testing. Traceability of all of the components used is assured and supporting documentation can be supplied upon request.

Quality Assurance Review Date: 3/9/11

Reviewed By Tracy M. Terry, QA Assistant: Tracy M. Terry

NOTE:

- Prior to receipt at Bionique Testing Laboratories, Inc., the stability of the test article is the responsibility of the company submitting the sample. Bionique Testing Laboratories Inc. will assume responsibility for sample stability following receipt and prior to being placed on test.
- This test is for the detection of microbiological growth and does not require statistical validation.

Page 1 of 2

APPENDIX

BIONIQUE TESTING LABORATORIES, INC.

Document ID #: DCF9002F
Title: QUALITY ASSURANCE REPORT - GMP
Effective Date: 03/12/10
Edition #: 01

REFERENCES

Regulatory:

- Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Code of Federal Regulations [CFR], Title 21 CFR Part 210, Current Good Manufacturing Practice in Manufacturing, Processing, Packing, or Holding of Drugs; General. FDA. Office of the Federal Register, National Archives and Records Department.
- Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Code of Federal Regulations [CFR], Title 21 CFR Part 211, Current Good Manufacturing Practice for Finished Pharmaceuticals. FDA. Office of the Federal Register, National Archives and Records Department.
- Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Points to Consider in the Characterization of Cell Lines Used to Produce Biologicals, Director, Center for Biologics Evaluation and Research, FDA. May, 1993. Docket No. 84N-0154.
- Department of Health and Human Services, Food and Drug Administration (USA) [FDA]. Code of Federal Regulations [CFR], Title 21 CFR Part 610.30, General Biological Products Standards; Subpart D, Test for Mycoplasma. FDA. Office of the Federal Register, National Archives and Records Department.

General:

- Barile MF, Kern J. Isolation of Mycoplasma arginini from commercial bovine sera and its implication in contaminated cell cultures. Proceedings of the Society for Experimental Biology and Medicine, Volume 138, Number 2, November 1971.
- Chen, T.R. In situ detection of mycoplasma contamination in cell cultures by fluorescent Hoechst 33258 stain. Experimental Cell Research, 104: 255-262, 1977.
- Carolyn K. Lincoln and Daniel J. Lundin. Mycoplasma Detection and Control. U. S. Fed. for Culture Collections Newsletter, Vol. 20, Number 4, 1990.
- Fetal Bovine Serum; Proposed Guideline. National Committee For Clinical Laboratory Standards (NCCLS), Vol. 10, Number 6, 1990. (NCCLS publication M25-P).
- McGarry GJ, Sarama J, Vanaman V. Cell Culture Techniques. ASM News, Vol. 51, No. 4, 1985.
- Tully JG, Razin S. Methods in Mycoplasma, Volumes I and II. Academic Press, N.Y., 1983.
- Barile MF, Razin S, Tully JG, Whitcomb RF. The Mycoplasmas, Volumes 1-4. Academic Press, N.Y., 1979.
- <http://www.bionique.com/> - Safe Cells Insights

Page 2 of 2

MYCOPLASMA TESTING SERVICES

APPENDIX IV

Document #: DCF3013D
Edition #: 10
Effective Date: 07/15/2003
Title: M-250 FINAL REPORT SHEET

Page 1 of 2

M-250 FINAL REPORT

Direct Specimen Culture
Procedures 3008, 3011, 3013

TO: WiCell QA
WiCell Research Institute
505 S. Rosa Rd., Suite 120
Madison, WI 53719
PHONE: 608-441-8019 FAX: 608-441-8011

BTL SAMPLE ID#: 64139 P.O.#: RP3891 DATE REC'D: 02/09/2011

TEST/CONTROL ARTICLE:

WA22-WB0056 #10027

LOT#: NA

DIRECT CULTURE SET-UP (DAY 0) DATE: 02/09/2011
INDICATOR CELL LINE (VERO) SEE DNA FLUORESCENCE RECORD SHEET

	DATE	
THIOGLYCOLLATE BROTH	DAY 7 + ⊕	02/16/2011
	DAY 28 + ⊕	03/09/2011
BROTH-FORTIFIED COMMERCIAL		
0.5 mL SAMPLE	DAY 7 + ⊕	02/16/2011
	DAY 28 + ⊕	03/09/2011
6.0 mL BROTH		
BROTH-MODIFIED HAYFLICK		
0.5 mL SAMPLE	DAY 7 + ⊕	02/16/2011
	DAY 28 + ⊕	03/09/2011
6.0 mL BROTH		
BROTH-HEART INFUSION		
0.5 mL SAMPLE	DAY 7 + ⊕	02/16/2011
	DAY 28 + ⊕	03/09/2011
6.0 mL BROTH		

(See Reverse)

APPENDIX IV

Page 2 of 2

Document #: DCF3013D
Edition #: 10
Effective Date: 07/15/2003
Title: M-250 FINAL REPORT SHEET

SAMPLE ID#	64139	AEROBIC	MICROAEROPHILIC	DATE
AGAR PLATES-FORTIFIED COMMERCIAL	DAY 7 + ⊕ DAY 14 + ⊕ DAY 21 + ⊕	+ ⊕ + ⊕ + ⊕	+ ⊕ + ⊕ + ⊕	02/16/2011 02/23/2011 03/02/2011
AGAR PLATES-MODIFIED HAYFLICK	DAY 7 + ⊕ DAY 14 + ⊕ DAY 21 + ⊕	+ ⊕ + ⊕ + ⊕	+ ⊕ + ⊕ + ⊕	02/16/2011 02/23/2011 03/02/2011
AGAR PLATES-HEART INFUSION	DAY 7 + ⊕ DAY 14 + ⊕ DAY 21 + ⊕	+ ⊕ + ⊕ + ⊕	+ ⊕ + ⊕ + ⊕	02/16/2011 02/23/2011 03/02/2011
BROTH SUBCULTURES (DAY 7)				DATE: 02/16/2011
AGAR PLATES-FORTIFIED COMMERCIAL	DAY 7 + ⊕ DAY 14 + ⊕ DAY 21 + ⊕	+ ⊕ + ⊕ + ⊕	+ ⊕ + ⊕ + ⊕	02/23/2011 03/02/2011 03/09/2011
AGAR PLATES-MODIFIED HAYFLICK	DAY 7 + ⊕ DAY 14 + ⊕ DAY 21 + ⊕	+ ⊕ + ⊕ + ⊕	+ ⊕ + ⊕ + ⊕	02/23/2011 03/02/2011 03/09/2011
AGAR PLATES-HEART INFUSION	DAY 7 + ⊕ DAY 14 + ⊕ DAY 21 + ⊕	+ ⊕ + ⊕ + ⊕	+ ⊕ + ⊕ + ⊕	02/23/2011 03/02/2011 03/09/2011

RESULTS: No detectable mycoplasma contamination

3/9/11
Date

Shayna E. Armstrong
Laboratory Director
Shayna E. Armstrong, Ph.D.

M-250 Procedural Summary: The objective of this test is to ascertain whether or not detectable mycoplasmas are present in an in vitro cell culture sample. As it is a primary culture, hybridoma, master seed stock or cell line. This procedure combines an indirect DNA staining approach to detect non-viable mycoplasmas with a direct culture methodology utilizing three different mycoplasma media formulations. The indirect approach involves the inoculation of the sample into a mycoplasma-free medium (DMEM) and performing a DNA fluorescence assay. The sample is inoculated into each of the 3 broth formulations and also onto duplicate plates (0.1 mL) for each of the 3 agar formulations. Subculture from broth to fresh agar plates is carried out after 7 days incubation. Agar plates are incubated aerobically and microaerophilically in order to detect any colony forming units morphologically indicative of mycoplasma contamination. Issues of the final report with signature of the Laboratory Director signifies that the required controls were performed concurrently with the test sample(s) as detailed in the referenced SOP and that all test conditions have been found to meet the required acceptance criteria for a valid test, including the appropriate results for the positive and negative controls.

Document ID #: DCF3008A
Title: **DNA FLUOROCHROME ASSAY RESULTS**
Effective Date: 3/24/10
Edition #: 07

DNA-FLUOROCHROME ASSAY RESULTS
Procedures 3008, 3009, 3011

Sample ID # **64139** **M-250** Date Rec'd: **02/09/2011** P.O. # **RP3891**

Indicator Cells Inoculated: Date/Initials: 2/15/11 / mk

Fixation: Date/Initials: 2/14/11 / mk

Staining: Date/Initials: 2/14/11 / mk

TEST/CONTROL ARTICLE:

WA22-WB0056 #10027

LOT# **NA**

WiCell QA
WiCell Research Institute
505 S. Rosa Rd., Suite 120
Madison, WI 53719

Phone: **608-441-8019**

Fax #: **608-441-8011**

DNA FLUOROCHROME ASSAY RESULTS:

☒ **NEGATIVE:** A reaction with staining limited to the nuclear region, which indicates no mycoplasmal contamination.

☐ **POSITIVE:** A significant amount of extranuclear staining which strongly suggests mycoplasmal contamination.

☐ **INCONCLUSIVE:**
☐ A significant amount of extranuclear staining consistent with low - level mycoplasmal contamination or nuclear degeneration.

☐ A significant amount of extranuclear staining consistent with bacterial, fungal or other microbial contaminant or viral CPE. Morphology not consistent for mycoplasmal contamination.

COMMENTS:

Date: 2/14/11 Results Read by: mk Date of Review: 2/14/11 Reviewed by: ll

Report Date: January 20, 2011

Case Details:

Cell Line: WA22-WB0056 10012

Passage #: 12

Date Completed: 1/20/2011

Cell Line Gender: Female

Investigator: Wisconsin International Stem Cell Bank

Specimen: hESC on Matrigel

Date of Sample: 1/12/2011

Tests, Reason for: Lot release testing

Results: 46,XX

Completed by Erik McIntire, CG(ASCP), on 1/20/2011

Reviewed and interpreted by Karen Dyer Montgomery, PhD, FACMG, on 1/20/2011

Interpretation: No clonal abnormalities were detected at the stated band level of resolution.



Cell: S02-43

Slide: 5(26)KARYOTYPE

Slide Type: Karyotyping

of Cells Counted: 40

of Cells Karyotyped: 4

of Cells Analyzed: 9

Band Level: 400-450

Results Transmitted by Fax / Email / Post

Sent By:

QC Review By:

Date:

Sent To:

Results Recorded: